


### **ashing muffle furnace with aluminum silicate inner chamber**

Ashing muffle furnace is a high-temperature experimental heating equipment, widely used in sample ashing, heat treatment, sintering, preheating, drying and other experimental operations, suitable for high-temperature calcination of inorganic and organic materials or residue analysis.

#### **Ashing Muffle Furnace**

The main features of the ashing muffle furnace:

1. high temperature capability: the working temperature can reach 1200°C, to meet the needs of most laboratory high temperature treatment.
2. integrated structure design: the whole machine is small in size, easy to place on the lab bench, saving space.
3. high-temperature resistant furnace chamber material: closed aluminum silicate furnace chamber, excellent heat preservation performance, fast warming, less heat loss.
4. Parallel left side opening structure: open the door from the left side parallel to the furnace door, to avoid the high temperature area close to the operator, improve safety.



5. intelligent temperature control system: the use of PID regulation technology, you can set the program through the key to keep the temperature stable.
6. perfect safety design: with over-temperature alarm, menu lock and open the door automatically power off and other multiple protection functions.
7. heating system optimization: nickel-chromium aluminum resistance wire embedded in the interior of the furnace chamber wall, to avoid direct exposure, to extend the service life.
8. auxiliary facilities: with observation window and exhaust gas exhaust port, easy to observe the experimental process and exhaust harmful gases.

### **Core Advantages**

1. high efficiency and energy saving: the use of aluminum silicate high-temperature insulation materials, low energy consumption, fast heating speed.
2. high safety: multiple automatic protection mechanism, operation more assured.
3. accurate temperature control: PID intelligent control technology, can realize constant temperature and warming curve control.
4. convenient operation: simple panel design and digital key input, easy to grasp.
5. adaptable: can be used in many types of laboratory and factory testing scenarios.

## Working Principle

Ashing Muffle Furnace is heated by converting electric energy into heat energy, using nickel-chromium-aluminum alloy heating wire, which generates heat and transfers it to the inner chamber after being energized. The furnace chamber is insulated with aluminum silicate high temperature insulation material to ensure that the heat is concentrated in the sample area. After the user sets the temperature through the PID controller, the system automatically adjusts the heating power to stabilize the furnace temperature at the set value. After the temperature is reached, the control system makes dynamic adjustments through the feedback loop to keep the heating curve smooth. In addition, the opening and closing of the furnace door is linked with the power supply to realize power failure protection and enhance the safety of use.

ashing muffle furnace with aluminum silicate inner chamber

| <b>Model</b>            | <b>AF201</b>                                | <b>AF202</b> | <b>AF203</b> |
|-------------------------|---|--------------|--------------|
| Door opening method     | Open the door parallel to the left          |              |              |
| Heating method          | heating on three sides, periodic heating up |              |              |
| Temperature range       | 300°C to 1200°C                             |              |              |
| Temperature resolution  | 1°C   |              |              |
| Temperature fluctuation | ±5%°C                                       |              |              |

ashing muffle furnace with aluminum silicate inner chamber

| <b>Model</b>                | <b>AF201</b>  | <b>AF202</b> | <b>AF203</b> |
|-----------------------------|---|--------------|--------------|
| Ambient to temperature time | 75 minutes  | 80 minutes   | 80 minutes   |
| Exhaust port                | diameter 22mmx1, rear chimney                                     |              |              |
| air inlet                   | diameter 15, middle of the door                                   |              |              |
| Inner Chamber material      | aluminum silicate   |              |              |
| Heating elements            | nickel chromium alloy wire, OCr27A17Mo2                           |              |              |
| Heating power               | 2.5Kw   | 5kw          | 10kw         |
| External structure          | cold rolled steel plate with chemically resistant surface coating |              |              |
| Temperature control method  | PID   |              |              |
| Program mode                | single segment program operation                                  |              |              |
| Temperature setting method  | Light touch four-button action with dual digital display          |              |              |
| Timer                       | 0 to 9999minutes  |              |              |

ashing muffle furnace with aluminum silicate inner chamber

| <b>Model</b>             | <b>AF201</b>   | <b>AF202</b>      | <b>AF203</b>      |
|--------------------------|--|-------------------|-------------------|
| Operating function       | Fixed value operation, program operation, timed operation, auto stop |                   |                   |
| Additional functionality | calibration function   |                   |                   |
| Safety device            | overcurrent protection switch  |                   |                   |
| Internal dimensions      | w120Xd200Xh80 mm   | w200Xd300Xh120 mm | w250Xd400Xh160 mm |
| External dimensions      | w420Xd658Xh610 mm  | w500Xd760Xh680 mm | w560Xd860Xh740 mm |
| Furnace chamber Volume   | 2L   | 7L                | 16L               |
| Power supply             | 230Vac, 50Hz, 11A  | 230Vac, 50Hz, 23A | 380Vac, 50Hz, 12A |
| Net weight               | 67kg   | 134kg             | 205kg             |
| Gross weight             | 80kg   | 150kg             | 230kg             |
| Standard configuration   | Chimney, open the door and cut off the electricity                   |                   |                   |